

SSAT 066: Statistical Analyses Results

Version 6.0: 22nd January 2018

eGFR - CKD-EPI with serum creatinine and cystatin C

mL/min/1.73 m ²	Total=60					
	Median (IQR)					
	Arm 1: <i>Truvada® (200mg/245mg film coated tablets) one tablet once daily + Raltegravir 400mg (one tablet) twice daily</i> N=20					
eGFR - CKD EPI Cr - Cys	92.0 (83.7 to 102.6) Range: 76.7 to 117.5	N=19 95.0 (86.2 to 102.8) Range: 78.5 to 115.7	N=19 97.9 (89.6 to 105.2) Range: 74.2 to 110.9	N=19 102.2 (97.3 to 110.7) Range: 67.5 to 123.0	N=19 111.8 (105.2 to 113.8) Range: 89.5 to 127.1	N=19 101.3 (88.6 to 107.4) Range: 67.9 to 118.5
	Arm 2: <i>Truvada® (200mg/245mg film coated tablets) one tablet once daily + Dolutegravir 50mg once daily</i> N=20					
eGFR - CKD EPI Cr - Cys	101.7 (93.8 to 107.2) Range: 63.1 to 121.9	88.9 (82.9 to 99.5) Range: 65.3 to 123.9	91.1 (85.1 to 95.8) Range: 56.0 to 114.6	95.0 (89.6 to 103.7) Range: 76.7 to 121.9	106.5 (97.2 to 110.7) Range: 90.3 to 115.7	100.7 (89.0 to 105.4) Range: 65.4 to 118.5
	Arm 3: <i>Stribild®150mg/150mg/200mg/245mg tablet once daily</i> N=20					
eGFR - CKD EPI Cr - Cys	109.1 (97.3 to 115.6) Range: 85.7 to 134.9	103.2 (91.4 to 113.0) Range: 77.6 to 133.0	104.1 (87.1 to 110.3) Range: 71.4 to 135.2	96.9 (90.0 to 112.3) Range: 74.8 to 146.6	111.4 (103.5 to 118.4) Range: 74.8 to 146.6	104.3 (94.6 to 114.8) Range: 70.0 to 140.3